

PLASTIC TIP SYSTEM

Abstract of Disclosure

A pre-packaged, sterilized, single-use dental scaler tip, to be supplied to dental professionals, is used with an insert for an ultrasonic dental scaler. The single-use tip is employed for maintenance and preventive care for patients with titanium and other dental implants. The tip is packaged individually and sterilized by gas or irradiation, only to be removed after purchase, such as by the dental professional just before use on a patient. Upon completion of the patient care, the tips are preferably discarded and the dental insert is autoclaved or otherwise sterilized.

Figures